

INTEGRAL ACTUATOR MODULES

FONCTIONS	E23041	E23042	E23045	E23046	E23171	E23271	E23150
ON-OFF switch						✓	
OUTPUT (motor)	1	2	1	2	1	1	4
DLO (Drive Lock Out)	✓	✓	✓	✓			✓
HS (Half Speed)		✓		✓			✓
HPWR (PTC)	4Amp	4Amp	6Amp	6Amp	6Amp	6Amp	6Amp
Programmable accel-decel	✓	✓	✓	✓	✓	✓	✓
Vspeed (programmable speed)							
EREAD Programmable current stop					✓	✓	✓
Stand-by mode					✓	✓	✓
Status LED (bicolor)						✓	✓
Attendant output	✓	✓	✓	✓	✓		✓