

PROGRAMMING KEYBOARD

E24001
E24002
E24003-SERIES

FUNCTIONAL DESCRIPTION

The programming keyboard includes a round CAN Bus connector (6 pins).

To use the programming keyboard, you only have to plug a wire to this CAN Bus connector and to a similar CAN Bus connector on any Liber-T Medtech powered up unit. Then turn ON the programming keyboard with the 0/1 switch.

SPECIFICATIONS

- Plug on LTM BUS
- General Settings
- 8 Pre-Programs
- History of Faults & Warnings
- Module Versions

