

Inflation... made easy!

i-PRESSURETM

• • • The perfect
air cushion companion



The i-PRESSURETM sensor is the ideal tool for EASY and TAILORED air cushion inflation. Using visual display and alarms, the sensor shows when to add or remove air from the cushion.

Compact design and robust housing, the i-PRESSURETM sensor features auto-on, extended tubing, cushion adaptability, disconnect and bottom-out alarm making it the perfect air cushion companion.

- Display Cushion Inflation
- Under & Over Pressure Alarms
- Auto-On
- Extended Tubing
- Disconnection Alarm
- Low Battery Alarm

www.i-PRESSURE.com


2855 De Celles, suite 300, Québec, QC G2C 1K7
T 1.418.842.2412 | F 1.418.845.6926 | info@i-pressure.com

LIBER • T
MEDTECH


i-PRESSURE™




SIMPLY FOLLOW THE LEDS



OVER INFLATION... REMOVE AIR



UNDER INFLATION... ADD AIR



INSTALLATION:



> HEAD LOOP



> TWO LOOPS



SPECIFICATIONS:

Lithium batteries duration	12 months (variable)
Battery type	2 AA
Housing	Black ABS plastic
Dimension	3.7" x 2.2" x 1.35"
Warning selection (audible or/and visual)	Slide switch inside battery compartment
Sound outputs (30 cm)	≥ 60 dB
Reverse polarity protection	Yes
Connect to	ROHO, STAR, others
Warranty*	1 year

* Refer to completed warranty for exclusions and limitations.