

USER INTERFACE

E30080 – B SWITCH

E30080 SERIES

FUNCTIONAL DESCRIPTION

The A-Switch series offers small, convenient remote control boxes to activate seat functions on a variety of rehab applications.

The **30080** model features a single toggle-type switch (momentary) that can be connected to Liber-T's E23040 single or dual integral actuator module and other product within Liber-T range.

Simple, costless and yet so efficient!



SPECIFICATIONS

- Black ABS plastic casing
- Waterproof toggle switch (MOM) – OFF – (MOM)
- Cable strain relief
- Jack audio plug 3.5mm straight
- Fixation: Metric or standard inserts
- Compatible TASH

DIMENSIONS

BLACK ABS PLASTIC CASING

Width	0.79"
Length	1.98"
Height	1.38"

