

USER INTERFACE

E30083 – DAS TOGGLE

E30083 SERIES

FUNCTIONAL DESCRIPTION

This series offers robust remote toggle switching to activate seat functions on a variety of rehab applications. A push button makes the ON/OFF functionality to output.

SPECIFICATIONS

- Black die cast aluminium casing
- Toggle switch (mom)-off-(mom) Hi-Amp
- Push-button switch on-off
- Jack audio plug 3.5mm connection (right angle)
- Fixation: Metric or standard inserts





DIMENSIONS

DIE CAST ALUMINIUM BLACK CASING	
Width	1.97"
Length	3.94"
Height (without switch stick)	0.90"
Height (with switch stick)	2.75"
Harness length	≅ 1750 mm

Jack audio connection (3.5mm)

SCHEMATICS



Positions

