

# USER INTERFACE

## E30090 – A-SWITCH

## E30090 SERIES

### FUNCTIONAL DESCRIPTION

The A-Switch series offers small, convenient remote control boxes to activate seat functions on a variety of rehab applications.

The **30090** switch box is used as a multi-function switch interface that features a 4-position joystick, 2 independent push buttons and a remote audio input.

Simple, costless and yet so efficient!



### SPECIFICATIONS

- Black ABS plastic casing
- 4 axis joystick switch
- 2 push switch button (red)
- 3.5 jack auxiliary (Forward Command)
- Cable strain relief
- DB-9 connector
- Fixation: Metric or standard inserts

### DIMENSIONS

#### BLACK ABS PLASTIC CASING

Width	2.13"
Length	3.26"
Height	1.20" – 2.70"

