

USER INTERFACE

E30095 – HEX SWITCH

E30090 SERIES

FUNCTIONAL DESCRIPTION

Those series offers small, convenient remote control boxes to activate seat functions on a variety of rehab applications.

The **E30095** switch box is used as a multi-function switch interface that features a 4-positions joystick. 2 LEDS (red and yellow) are visible to advise the user. If you are connected to the Quattro (E23050) the DLO will be red and the HS, yellow.

SPECIFICATIONS

- Black ABS plastic casing
- 4 axis joystick switch
- 2 LEDS (red, yellow)
- DB9 female connection
- Fixation: Metric or standard inserts
- Customized faceplate available



DIMENSIONS

PLASTIC ABS BLACK CASING

Width	2.13"
Length	3.26"
Height (without switchstick)	1.20"
Height (with switchstick)	3.20"
Harness length	≅ 5.00'

DB9 CONNECTION

# PIN	FUNCTION
1	Forward joystick
2	Reverse joystick
3	Right joystick
4	Left joystick
5	Yellow LED (cathode)
6	Red LED (cathode)
7	Not connected
8	Common joystick
9	GND LEDs (commun anodes)

